

Modern Physics

For the complete bibliography with references as well as an explanation of the classification scheme go to:

[Demonstration Bibliography](#)

The **demonstration name** listed in the bibliography is either the name listed on the reference or, if none is given, a simple descriptive name. In cases where there are several common names for a demonstration, the committee has chosen a preferred name.

The **description** is very brief. It is not intended to be a summary of the reference. One sentence is, in general, sufficient to describe the unique characteristics, if any, of an item. Each source has a unique numbering format. These unique formats are used identify references in the Bibliography.


The formats for the **reference** column and links to the sources are listed below:

Reference	Source
M-1	Sutton
Ma-1	Freier & Anderson
M-1d	Hilton
8-2.8	Meiners
M-108	Dick & Rae
1A 12.01	University of Minnesota Handbook
AJP 52(1),85	American Journal of Physics
TPT 15(5),300	The Physics Teacher
Disc 01-01	The Video Encyclopedia of Physics Demonstrations
PIRA 200	Physics Instructional Resource Association
PIRA 500	PIRA 500
PIRA 1000	PIRA 1000


Each demonstration is listed in only one location, even if it is commonly used to illustrate several concepts.

KINETIC THEORY	4
Brownian Motion	4
Brownian Motion Cell	4
QUANTUM EFFECTS	5
Photoelectric Effect.....	5
Discharging Zinc Plate.....	5
Stopping Potential.....	6
Millikan Oil Drop	7
Millikan Oil Drop Experiment.....	7
Wave Mechanics	8
Frustrated Total Internal Reflection.....	8
Vibrating Circular Wire	9
X-ray and Electron Diffraction	10
Electron Diffraction	10
Microwave Bragg Diffraction.....	11
Sample X-Ray Tube.....	12
Condensed Matter	13
Electron Conduction Model.....	13
PHOTOMETRY	14
Blackbodies.....	14
Hole in a Box	14
Black Body Spectrum	15
ATOMIC PHYSICS	16
Spectra.....	16
Student Gratings and Line Sources.....	16
Absorption.....	17
Band Absorption Spectrum.....	17
Resonance Radiation.....	18
Projected Mercury Spectrum	18
Fluorescence and Phosphorescence	19
Particle in a Box	20
Ionization Potentials.....	21
Plasma Globe	21
Atomic Models.....	22
Electron Orbital Models.....	22
NUCLEAR PHYSICS	23
Radioactivity.....	23
Geiger Counter and Samples	23
Half Life with Isotope Generator.....	24
Half Life of Silver.....	25
Cosmic Rays	26
Conduction in Gases	27
Discharge an Electroscope.....	27
Nuclear Reactions	28
Mousetrap Chain Reaction.....	28
Particle Detectors	29

Ludlum Detectors.....	29
Deflection of Beta Rays	30
Diffusion Cloud Chamber.....	31
Models of the Nucleus	32
Rutherford Scattering.....	32
Rutherford Scattering on an Air Table.....	33
Rutherford Model.....	34
RELATIVITY	35
Special Relativity	35
Induction Coil Relativity.....	35
Lorentz Transformation / Time Dilation.....	36

Brownian Motion**Brownian Motion Cell**


Observe the motion of particles in a smoke cell through a microscope. The balls are about 1 micron in diameter. Put a drop of the solution on a slide then place a cover slide on top of it. Use a drop of oil with the oil immersion objective. The well slide works best (less systematic movement); if you use it, place a drop in the well then dilute it with water.




Location: Ga3

Photoelectric Effect

Discharging Zinc Plate




A polished (with fine abrasive paper) zinc plate sits atop an electroscope. Charge the electroscope with a rod. Turn on the UV light. Expose the zinc plate to the UV rays. Do not turn the UV light towards the students. Compare what happens when you charge the plate positively with what happens when you charge it negatively.



Location: Ka2, Gc2

MODERN PHYSICS	7A10.30	QUANTUM EFFECTS
Photoelectric Effect		
<h1>Stopping Potential</h1>		
	h/e apparatus. Measure the stopping potential of the lines of the mercury spectrum with a photodiode.	



Location: Ka2

Millikan Oil Drop**Millikan Oil Drop Experiment**

The oil drop experiment.




Location: Ka3

Wave Mechanics


Frustrated Total Internal Reflection

Squeeze two right angle prisms together using a C Clamp and some notched blocks of wood while directing a beam of light at the interface. Use the blackboard optics equipment for this.




Location: Blackboard Optics Cart

Vibrating Circular Wire




The PASCO signal generator drives the PASCO mechanical vibrator. The three node resonance is at 20Hz. Turn up the amplitude after finding the resonance. Other resonances are at 60, 100, and 140 Hz. Pin the wire at 180 degrees from the mechanical vibrator to get full numbers of waves.




Location: Ka3, Ec2, Fb4

X-ray and Electron Diffraction

Electron Diffraction



Rings or spots diffraction patterns can be shown to the class using a TV camera. The beam and bias will already be adjusted. Turn up the HV until the pattern appears. Do not let the current exceed 170 microamps. Turn down the HV when you are done. The small jack at the back has a negative bias (20 V or so) on it. The filament voltage is 6.3 V and one side is connected to the common ground for the HV and negative bias supply.




Location: Ka4

X-ray and Electron Diffraction

Microwave Bragg Diffraction

Microwave diffraction is observed from a crystal model made of steel ball bearings mounted in a Styrofoam cube.




Location: Ia7, Ka5

X-ray and Electron Diffraction

Sample X-Ray Tube

Show a large X-ray tube.



Location: Kd1


Condensed Matter**Electron Conduction Model**

The Air Cushion Table is used to show bound charge carriers in an insulator; a free charge in an insulator, conduction in a semiconductor (N type and P type. Air Cushion Table manual: 2.4.9, 2.4.10, 2.4.11, 2.4.12, 2.4.13.




Location: Ga5

OPTICS	6B40.20	PHOTOMETRY
Blackbodies		
Hole in a Box		
<p>A box that is painted black has a hole in a hinged door. The hole appears blacker than the box, When the door is opened, you see the box is actually painted white on its inside surfaces.</p>		



Location: Jb2


OPTICS	6B40.40	PHOTOMETRY
Blackbodies		
<h1>Black Body Spectrum</h1>		
<p>A spectrometer and bolometer connected to Science Workshop is used to record the spectrum from an incandescent bulb at various temperatures.</p>		




Location: Science Workshop Cabinet, Optics Table Cabinet

Spectra

Student Gratings and Line Sources



Pass out the 1"x1" gratings to the students. These have 13,400 lines per inch. Turn on one of the light sources. There is a single filament white light source and three discharge tubes, Hg, He, and Ne. The Didymium filter can be placed in front of the white light source to show selective absorption. Also shown is a spectrometer with a TV camera showing hydrogen Balmer lines and the Ride Tide USB spectrometer showing the line spectrum on the computer. Use the door lite curtains to darken the room.




Location: Ja1, Jb2

Absorption

Band Absorption Spectrum

Didymium glass is used to show band absorption. The Reveal 60 bulb is used with both a grating and the Red Tide spectrometer to show the bands. An ordinary bulb is shown for comparison.




Location: Jb2, Ib1

Resonance Radiation

Projected Mercury Spectrum

The UV lines of the projected spectrum are made visible using a fluorescent board



Location: Ka2

Resonance Radiation

Fluorescence and Phosphorescence

Fluorescent and phosphorescent materials are shown with a black light.



Location: Kb3

Resonance Radiation

Particle in a Box

Nanobeads (quantum dots) of 4 different sizes in solution fluoresce when illuminated with a UV light. The energy calculation is based on the calculation of a particle in a potential sphere (see the appendix of the instructions for more info). The quantum "sphere" or "box" is not empty but filled with a semiconductor. In order for the students to understand this they must first start with the derivation of a 1D particle in a box calculation, the sphere is basically the same physical concept with a different mathematical geometry.

Using the Red Tide spectrometer the peaks from the different fluorescing vials can be measured and compared with theory.



Location: Kb3

Ionization Potentials

Plasma Globe

Commercial plasma tubes are discussed. Bring your hand near the globe.




Location: Kb4

Atomic Models

Electron Orbital Models



A set of Klinger electron orbital models.




Location: Kb5

Radioactivity

Geiger Counter and Samples



Turn on the counter and place sample under Geiger Tube. Test to see how distance from the sample affects the count. Place objects like: paper, wood, lead and other items between the source and the tube to check their shielding properties. The alpha paddle can be used to selectively detect alpha particles only.




Location: Kb6

Radioactivity

Half Life with Isotope Generator

Follow the instructions on how to “milk” the generator. It is good practice to survey the work area before and after for contamination. The Geiger counter can be connected to the computer and using Science Workshop one can plot the count rate as a function of time.



Location: Kb6

MODERN
PHYSICS

7D10.33

NUCLEAR PHYSICS

Radioactivity

Half Life of Silver

Measure the half life of silver activated by a neutron source. Use a glove to handle the “hot” silver as it is a beta emitter. Leave as the closest position to the neutron source for a few minutes to activate it.



Location: Advanced Lab, Kb6

MODERN PHYSICS	7D10.80	NUCLEAR PHYSICS
Radioactivity		
Cosmic Rays		
	Scintillator detectors show the passage of a cosmic ray by coincidence counting.	

Electricity and Magnetism	5D40.30	RESISTANCE
Conduction in Gases		
Discharge an Electroscope		
<p>Charge the electroscope either positively or negatively by induction (this gives you control over the charge which is necessary as the leaf electroscope is sensitive and can be easily overcharged). Then place the Po-210 source (0.5 mCi when new) in the electroscope and watch the electroscope discharge in a couple of seconds. Use safety protocols when handling the source.</p>		



Location: Kb6, Gc2

Nuclear Reactions**Mousetrap Chain Reaction**


A 72 mousetraps with corks set on them. Trigger with a single “neutron” (cork). This requires about 20 minutes and nerves of steel to set up. Each trap stores about a joule of energy when set.



Location: KbT

Particle Detectors**Ludlum Detectors**

Ludlum hand held alpha, beta and gamma detectors are used with a variety of sources.



Location: Kb6

MODERN PHYSICS	7D30.30	NUCLEAR PHYSICS
Particle Detectors		
Deflection of Beta Rays		
<p>Turn on the counter and lay it on its side. Open the beta window. Place a beta source nearby, being sure the emitting side is toward the counter. Then bring a strong magnet near the source. The count rate will drop.</p>		



Location: Kb6, Hd5

MODERN
PHYSICS

7D30.60

NUCLEAR PHYSICS

Particle Detectors

Diffusion Cloud Chamber

This demo requires advance notice to obtain ice and let the chamber cool (about 40 minutes). Tracks are most impressive if the students are allowed to view individually. The TV works quite well too. Use either ambient cosmic rays or Pb 210 as a source.



Location: Kc1

Models of the Nucleus

Rutherford Scattering

Balls roll down a ramp onto a potential surface to model Rutherford scattering.



Location: Kc2

Models of the Nucleus

Rutherford Scattering on an Air Table

The Air Cushion Table is set up with a plastic platform that holds a magnet which repels a magnetic puck that is set in motion beneath it (Air Cushion Manual, 2.5.1). One can also illustrate the scattering of alpha particles through a foil using a foil made of a 1D lattice four magnets that repel a magnetic puck that is set in motion toward it (Air Cushion Table Manual 2.5.2).



Location: Ga4

Models of the Nucleus

Rutherford Model

The Air Cushion Table is set up with a plastic platform that holds a magnet which attracts a magnetic puck that is set in motion beneath it (Air Cushion Manual, 2.5.3). The moving puck goes into orbit beneath the attracting magnet.




Location: Ga4

Special Relativity

Induction Coil Relativity

On using the simple induction coil and galvanometer as a special relativity demonstration: AJP 48(9), 780.



Location: BdT, Hc2, Ib3

Special Relativity

Lorentz Transformation / Time Dilation



The Mechanical Universe, chapter 42, and the Hewitt film "Relativistic Time Dilation"



Location: Ad5